

Carqon

# Cushion extension straps

Hereby you will receive the cushion extension-straps for your first Carqon Cruise bikes. With these straps you will be able to connect the seat-cushion to the backside of the box in a proper way.

**Please note;** for future orders, these straps will be automatically delivered together with your Carqon Cruise. In order to make the current accessory line for the Carqon Classic as much as possible compatible with the Carqon Cruise, we decided to come up with a solution to make the two-seater seat cushion fit well for both bikes.

- 1 Pull the original strap of the seat-cushion through the extend strap. Make sure the wide part of the strap is facing the cushion, and the silver part of the button is facing upwards.



- 2 Click the buttons together so that the two straps are fastened together.



**More information?**  
go to [support.carqon.com](https://support.carqon.com)



- 3 The left strap of the cushion is longer, which you can pull through the extend strap so the Carqon logo is visible again.



- 4 Now both straps are long enough in order to mount them properly on the box of the Carqon Cruise.



Enjoy your ride!

In case of questions, please contact our dealer support by e-mail at [info@carqon.com](mailto:info@carqon.com) or by phone +31 (0)33 74 10 743.

Best regards,

Team Carqon



**More information?**  
go to [support.carqon.com](https://support.carqon.com)

